Lastodur®

Strong/Soft

Langzugbinde
Long-stretch bandage
Bande à allongement long
Zwachtel met lange rek
Benda a lunga trazione
Venda de tracción larga
Ligadura de longa tração
Ελαστικός επίδεσμος υψηλής έκτασης
Dlouhotázně obinadlo
Dlhoťažné ovínadlo
Bandaż o dużej rozciągłości
Hosszú megnyúlású pólya
Эластичный бинт
Dolgo-elastíchni povoj
Långsträcksbinda
Langstræksbind
Pitkävennytteinen side
Langstrekksbind
Uzun germe bandaj
長拉力繃帶
رباط قابل للمط الى طول بالغ

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EN Instructions for use

General description of the device

The long-stretch compression bandages are narrow woven. They are provided with an end fixation. The bandages themselves are made of natural and synthetic materials in various proportions and are brown in colour. They don't have a cohesive layer. The non-active and non-sterile devices are available in different widths and lengths. Lastodur straff/strong and weich/soft are permanently elastic long-stretch bandages.

Properties

Extensibility: approx. 160 %.
Washable up to 15 times: Wash the bandages in regular wash cycle at a temperature of 60 °C and dry them flat. Do not bleach, tumble dry or dry clean the bandages.

Intended purpose and indications

The compression bandages are non-sterile medical devices used by healthcare professionals on intact skin of humans (above 6 months of age) in clinical and home environment. The compression bandages can be used for compression therapy (ABPI 0.9-1.3) (decongestion phase) for the treatment of:

- varicose veins (e.g. varicose in pregnancy, supportive in the invasive treatment of varicoses)
- chronic venous insufficiency (CVI) with or without venous leg ulcers (VLU) (e.g. CVI according to CEAP classification C3-C6, primary and secondary prevention of leg ulcer venous, venous insufficiency and angiodysplasia)
- thrombosis/thrombophlebitis (z. B. thrombophlebitis as well as condition after healed thrombophlebitis, deep leg and arm vein thrombosis, condition after deep vein thrombosis, post-thrombotic syndrome, thrombophlebitis)
- chronic oedema of other origin (e.g. oedema in pregnancy, post-traumatic oedema, hormonal oedema, lipodema, stasis conditions due to immobility (arthrogenic congestion syndrome, paresis), medically induced oedema [e.g. calcium antagonists, isosorbide dinitrate, lithium ointment, sex hormones])
- lymph oedema

The devices can also be used as:

- support and relief bandages in musculoskeletal injuries
- support and relief bandages in the management of limbs following amputation
- support bandages for distortions, contusions, tendon injuries, and dislocations (joints)

The compression bandages can also be used in combination with primary and secondary dressings. Lastodur straff/strong and weich/soft are reusable devices.

Composition

Lastodur straff/strong: 84 % cotton, 6 % polyamide, 10 % elastane

Lastodur weich/soft: 82 % cotton, 13 % polyamide, 5 % elastane

Kontraindications

The long-stretch compression bandages must not be used in case of:

- advanced peripheral arterial occlusive disease (PAOD)
- decompensated heart failure (NYHA III + IV)
- phlegmasia coerulea dolens
- ankle-brachial pressure index (ABPI) > 1.3 and ≤ 0.5
- septic phlebitis
- florid erysipelas

Precautions

Special precautionary measures are required in the event of:

- pronounced sensibility disorders of the extremities (e.g. in the case of advanced diabetic peripheral polyneuropathy)
- compensated peripheral arterial occlusive disease (PAOD) with a reduced ankle-brachial pressure index (ABPI 0.6–0.8)
- chronic compensated cardiac insufficiency

• incompatibility with bandage material

The application technique, e.g. with regards to stretch and/or padding, may have to be adapted to the patient's anatomical and pathophysiological characteristics, especially in case of fragile or aged skin.

If the foot pulse is not palpable, a Doppler ultrasound should be carried out to determine the ABPI prior to starting therapy. The right measurement method must be considered for the ABPI. During the use of the compression bandages, the ABPI must be controlled regularly. The therapy must be adapted accordingly.

The devices should not be used on infants below the age of 6 months. When applying the compression bandages on children (from 6 months on), the pressure should be reduced appropriately. If at any time during therapy the patient develops severe pain, numbness, tingling or noticeable discolouration of the toes, a healthcare professional must be contacted immediately, and the compression bandages must be removed.

Further remarks

- Due to high resting pressure, the bandages should be removed overnight.
- Do not let the bandages get wet during showering/bathing.
- Do not use the bandages as a primary dressing or directly on open wounds.
- Bandage clips are only used to fasten the end of the bandage for storage. Not to be used on patients.
- A slipping of the bandage indicates a loss of compression. In this case, a medically trained person has to reapply the bandage.

Incident reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/745 on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Product disposal

In order to minimize the risk of potential infection hazards, or environmental pollution, disposable components of the medical device should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards. Dispose of the medical device with regular hospital waste.

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