

UROMED Füllmedium

UROMED Filling Medium

REF 1550

Bringen Sie Konstanz in die Ballons Ihrer Silikonkatheter!

Create constancy in the balloons of your silicone catheters!

Füllmedium zum Füllen von urologischen Ballonkathetern

- enthält 10 % Glycerin in aqua dest.
- hält das Füllvolumen konstant
- erleichtert die Katheteranlage
- latexfrei

Filling medium for filling urological balloon catheters

- contains 10% glycerin in aqua dest.
- keeps the filling volume constant
- easier application of catheters
- free of latex



Diffusion bei Silikonballonkathetern

Die Ballons der Silikonballonkatheter lassen allgemein materialbedingt eine Diffusion zu. Dadurch kann sich der Ballon langsam entleeren. Das gilt im Besonderen für Nephrostomie- und suprapubische Ballonkatheter mit geringem Füllvolumen. Das UROMED Füllmedium hilft durch das enthaltende Glycerin eine Diffusion effektiv zu reduzieren.

Diffusion of silicone balloon catheters

The balloons of silicone balloon catheters in general allow a diffusion caused by the material. Therefore the balloon might evacuate slowly. This applies especially to nephrostomy and suprapubic balloon catheters with small filling volume. The UROMED filling medium helps to reduce any diffusion effectively due to the glycerin.

REF REF	Volumen Volume	PZN CPN
1550.05	05 ml	00730833
1550.10	10 ml	00730827

Hilfsmittelkatalog: 15.99.99.0006

Verpackung: steril
Originalkarton: 10 Stück
Packaging: sterile
Original carton: 10 pieces