



Product description

LoFric® Hydro-Kit™, single use, hydrophilic urinary catheter with salt solution for activation, intended for intermittent urinary catheterization.

Sterilization and shelf life

Sterilized by irradiation. Use before expiry date on package.

Package

30 catheters per customer box. 2 customer boxes per shipping box.

Material composition

Hydrophilic surface coating: Polyvinyl pyrrolidone (PVP).

Catheter: Polyolefin-based elastomer (POBE).

Sachet: Laminate foil of polyethylene terephthalate (PET), aluminum and polyethylene (PE). Sterile water with sodium chloride (NaCl).

Item packaging: Polyethylene butylacrylate (EBA).

Product range

Reference	Type	Length	Size
421 08 XXX*	LoFric® Hydro-Kit™ Nelaton (paediatric)	20 cm	CH 08
421 10 XXX*	LoFric® Hydro-Kit™ Nelaton (paediatric)	20 cm	CH 10
423 08 XXX*	LoFric® Hydro-Kit™ Nelaton (female)	20 cm	CH 08
423 10 XXX*	LoFric® Hydro-Kit™ Nelaton (female)	20 cm	CH 10
423 12 XXX*	LoFric® Hydro-Kit™ Nelaton (female)	20 cm	CH 12
423 14 XXX*	LoFric® Hydro-Kit™ Nelaton (female)	20 cm	CH 14
423 16 XXX*	LoFric® Hydro-Kit™ Nelaton (female)	20 cm	CH 16
423 18 XXX*	LoFric® Hydro-Kit™ Nelaton (female)	20 cm	CH 18
420 08 XXX*	LoFric® Hydro-Kit™ Nelaton (male)	40 cm	CH 08
420 10 XXX*	LoFric® Hydro-Kit™ Nelaton (male)	40 cm	CH 10
420 12 XXX*	LoFric® Hydro-Kit™ Nelaton (male)	40 cm	CH 12
420 14 XXX*	LoFric® Hydro-Kit™ Nelaton (male)	40 cm	CH 14
420 16 XXX*	LoFric® Hydro-Kit™ Nelaton (male)	40 cm	CH 16
420 18 XXX*	LoFric® Hydro-Kit™ Nelaton (male)	40 cm	CH 18
425 10 XXX*	LoFric® Hydro-Kit™ Tiemann (male)	40 cm	CH 10
425 12 XXX*	LoFric® Hydro-Kit™ Tiemann (male)	40 cm	CH 12
425 14 XXX*	LoFric® Hydro-Kit™ Tiemann (male)	40 cm	CH 14
425 16 XXX*	LoFric® Hydro-Kit™ Tiemann (male)	40 cm	CH 16
425 18 XXX*	LoFric® Hydro-Kit™ Tiemann (male)	40 cm	CH 18

*Generic article number without the three digit suffix specific for a region or country destination when distributing an article.



User-friendly packaging.



LoFric® Hydro-Kit™

References

1. Waller et al., Spinal Cord, 1997
2. Dunger et al., Acta Endocrinologica, 1988
3. Lundgren et al., Spinal Cord, 2000
4. Newman et al., Springer, 2018

At Wellspect we develop innovative continence care solutions that improve quality of life for people with bladder and bowel problems. We inspire our users to build self-confidence and independence as well as good health and well-being. We have been leading the industry for over 40 years with our product brands LoFric® and Navina™. We always aim to minimize the environmental impact of our products and passionately strive to become climate neutral. We work together with users and healthcare professionals to improve clinical outcome in a sustainable way, now and for the future.

Wellspect. A Real Difference.

For more information about our products and services, please visit [Wellspect.com](https://www.wellspect.com).

Join the conversation on LinkedIn, Twitter, Facebook and Instagram.

[wellspect.com](https://www.wellspect.com)

Wellspect HealthCare, Aminogatan 1, P.O. Box 14, SE-431 21 Mölndal, Sweden.
Phone: +46 31 376 40 00. www.wellspect.com

LoFric



Manufacturer: Wellspect HealthCare, Aminogatan 1, P.O. Box 14, SE-431 21 Mölndal, Sweden.
Phone: +46 31 376 40 00. www.wellspect.com

CE 2797

Rx Only